
<b>DIGITAL MICROSCOPE</b>  User Manual 330H	<b>CONTENTS</b>	<b>PREFACE</b> There are some instructions in this device. This product is complete to use, so please read the manual in detail to understand how to operate the microscope before use.	<b>DETAILLED EXPLANATION OF THE BUTTONS AND FUNCTIONS</b>	<b>PRODUCT DESCRIPTION</b>	<b>MODEL INTRODUCTION</b>	<b>CHARGING INSTRUCTIONS</b>	<b>TECHNICAL PARAMETERS</b>	<b>TROUBLESHOOTING</b>
	English.....01-10 Deutsch.....11-20 Italiano.....21-30 Français.....31-40 Español.....41-50 日本語.....51-60 中 文.....61-70	1. Please fully charge the device before use. Do not charge directly via PC. 2. Please open the cover before use. Close the cover after use. 3. The optimal focal length of the microscope is 110mm. You need to rotate the objective lenses to get the best view. 4. Do not disassemble the product or its components to avoid damage of the product.  <b>Product Description</b> Mode Switch Button (Mode Switch/Zoom) In observation mode, short press the Mode switch button to switch between: - Infrared, short press for infrared. - TV Out, short press for TV Out. - Microscope Stand.  OK Button (Image Selection/Confirm) In observation mode, short press the OK button to take a photo. TV Out (Video Output) In observation mode, long press to turn on. Microscope Stand (Microscope Stand) In observation mode, short press for the video setting page, and short press again to end recording. Settings Button (Setting/Delete) In observation mode, short press to enter the photo settings page, and short press again to exit. TV Out (Video Output) In observation mode, short press to turn on. Microscope Stand (Microscope Stand) In observation mode, short press for the video setting page, and short press again to end recording. OK Button (Image Selection/Confirm) In observation mode, short press the OK button to take a photo. TV Out (Video Output) In observation mode, short press to turn on. Microscope Stand (Microscope Stand) In observation mode, short press for the video setting page, and short press again to end recording. Settings Button (Setting/Delete) In observation mode, short press to enter the photo settings page, and short press again to exit. OK Button (Image Selection/Confirm) In observation mode, short press the OK button to take a photo. TV Out (Video Output) In observation mode, short press to turn on. Microscope Stand (Microscope Stand) In observation mode, short press for the video setting page, and short press again to end recording. Settings Button (Setting/Delete)	1. <b>OK Button (Image Selection/Confirm)</b> In observation mode, short press the OK button to take a photo. 2. <b>Power (TV Out)</b> In observation mode, long press to turn on. 3. <b>Microscope Stand</b> In observation mode, short press to turn off the objective lenses. 4. <b>Mode Switch Button (Mode Switch/Zoom)</b> In observation mode, short press the Mode switch button to switch between: - Infrared, short press for infrared. - TV Out, short press for TV Out. - Microscope Stand.  <b>Product Description</b> Mode Switch Button (Mode Switch/Zoom) In observation mode, short press the Mode switch button to switch between: - Infrared, short press for infrared. - TV Out, short press for TV Out. - Microscope Stand.  OK Button (Image Selection/Confirm) In observation mode, short press the OK button to take a photo. TV Out (Video Output) In observation mode, short press to turn on. Microscope Stand (Microscope Stand) In observation mode, short press for the video setting page, and short press again to end recording. Settings Button (Setting/Delete) In observation mode, short press to enter the photo settings page, and short press again to exit. TV Out (Video Output) In observation mode, short press to turn on. Microscope Stand (Microscope Stand) In observation mode, short press for the video setting page, and short press again to end recording. OK Button (Image Selection/Confirm) In observation mode, short press the OK button to take a photo. TV Out (Video Output) In observation mode, short press to turn on. Microscope Stand (Microscope Stand) In observation mode, short press for the video setting page, and short press again to end recording. Settings Button (Setting/Delete)	1. <b>Windows Camera</b> Windows camera, as shown below.	Opening the software recommended by the "Photo Booth" software with your Apple computer:	1) Firstly, connect the microscope and the screen via an HDMI cable to the computer, and then turn on the power of the microscope. Connect the microscope to the computer for card reading writing, and then charge the battery. Turn on the power of the microscope, and then charge the battery. Turn on the power of the microscope, and then charge the battery.	Screen parameters * IPS Pixels 5.0 megapixels Magnification 1200X Image format JPG Video format AVI Battery capacity 2000 mAh Operating time 3 hours Charging time 1 hour Supported system Windows vista/XP/7/10/11 MacOS X 11.0 or higher	Q1: The device cannot be charged. A1: Please charge the device via the power supply. Q2: Turn on the main power of the power supply (e.g., power strip, computer power). Q3: Turn on the power of the main unit with a line needs to force shutdown the device.

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